00203

Allegato 6 della I71201/PT rev. 0

Linee guida
PT 2026-Fitofarmaci e Nitrati
ARPAE

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1. FOREWORD

The purpose of the PT is to verify the performance of participants in accordance with the current revision of the SANTE *"Guidance document on analytical quality control and method validation procedures for pesticides residues and analysis in food and feed"* which provides for an extended uncertainty of 50% for pesticides. For the evaluation of the performances of nitrates a standard deviation of 12% is used equal to an extended uncertainty of 24% deduced from inter- comparison collected over time.

ARPAE has been accredited as a proficiency test provider since 25/09/2019, in compliance with the UNI CEI EN ISO/IEC 17043:2024 standard, to which the PTP code N°: 00203P has been assigned.

This guideline, and all the documentation sent to the participants, are an integral part of the QMS documents, drawn up in compliance with the UNI CEI EN ISO/IEC 17043:2024 standard.

The proficiency tests are organized by ARPAE Ferrara for pesticide residues on fruit and vegetables and wine, with the application of multi-residue analysis methods, and for nitrates on leafy vegetables. The organization of the proficiency tests aims to improve the quality, accuracy and comparability of data on pesticide residues and nitrate content; as well as to evaluate the ability of each laboratory to operate with results that cover the entire range of maximum residues allowed, in all groups of fruit and vegetable matrices.

For the current year, ARPAE is planning to organize 3 proficiency sessions tests , each of which includes both the preparation of samples for the determination of pesticide residues and samples for the determination of nitrates.

All analytical determinations concerning the test samples are carried out by an accredited laboratory according to the Arpae UNI CEI EN ISO/IEC 17025:2018 standard .

2. CONFIDENTIALITY

The laboratories are identified only and exclusively through a numerical code.

The codes are assigned on a random basis and communicated to the participant when the Google Sheet is shared for the transmission of their results.

In the event that the results of a participant be communicated to third parties, this occurs only with the prior knowledge and written authorization of the participant.

This authorization is not foreseen in cases where the request is required by law, however the Provider is required to give written notice to the participant.

	<p>Linea guida Fitofarmaci e Nitrati Pesticides and Nitrates Guideline</p>	<p>Allegato 6 I71201/PT rev.0</p>
---	--	---------------------------------------

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3. **PESTICIDES TEST SAMPLES**

Preparation for the proficiency test, (PT), includes:

- Choice of matrix from those present in the table in Annex A of the SANTE "*Guidance document on analytical quality control and method validation procedures for pesticides residues and analysis in food and feed*". Annually, over the three sessions, depending on the availability of the matrices, we will try to choose a vegetable product, a fruit product and one with a wine matrix.
- Choice of active substances for the preparation of the boosting solution. The choice is made by comparing the active substances present in Annex 4 with the RASFF "Rapid Alert System for Food and Feed ", the DIGISAN operational guidelines "General Directorate for Hygiene, Food Safety and Nutrition" of the Ministry of Health, and the fields of use on the chosen matrix.
- Homogenization of an adequate amount of matrix.
- Mixture with spiked active substances.
- Subdivision into appropriate containers.

Ten such bins, chosen at random, are analyzed to evaluate homogeneity.

Samples are kept frozen before shipping to participants and shipped on dry ice.

Another six containers, again chosen at random, are subsequently analyzed in duplicate, two at the time of shipping of the samples, two 48 hours after shipping of the objects and with the portion kept at refrigerated temperature and two at the time of the dead-line, for evaluate stability.

4. **NITRATE TEST SAMPLES**

The preparation of the PT involves the homogenization of an adequate quantity of matrix and the division into appropriate containers. Ten such containers, chosen at random, are analyzed to evaluate homogeneity.

Samples are kept frozen before shipping to participants and shipped on dry ice.

Another six containers, again chosen at random, are subsequently analyzed in duplicate, two at the time of shipping of the samples, two 48hours after shipping of the objects and with the portion kept at refrigerated temperature and two at the time of the dead-line, for evaluate stability.

In the absence of analytical data of homogeneity and/or stability, as required by point 6.1.1 of the ISO 13528:2022 Standard, the PTP can perform an evaluation based on data from previous similar rounds. The use of this approach will be highlighted in the Final Report(s)

	<p>Linea guida Fitofarmaci e Nitrati Pesticides and Nitrates Guideline</p>	<p>Allegato 6 I71201/PT rev.0</p>
---	--	---------------------------------------

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5. **COMMUNICATIONS**

Communications between ARPAE and PT participants take place using the following Annexes:

- Annex 1: PT Program for 2026 - sent by ARPAE via e-mail with the presentation letter.
- Annex 2: COSTS Proficiency Test 2026 - sent by ARPAE via email with the presentation letter.
- Annex 3: PT membership form for 2026 - sent by ARPAE via e-mail with the presentation letter, it must be returned completed within the indicated times
- Annex 4: list of active substances.
- Google Sheet: is prepared for each individual participant. For each session, information relating to the time taken to deliver the samples and the status upon arrival must be entered and the results and analytical methods adopted must be entered.

If during the year there are changes relating to company names, shipping addresses, billing details, email addresses and reference personnel, it is the responsibility of the participating laboratory to promptly notify the coordinator by completing a new Annex 3: PT membership card year 2026.

6. **PRODUCTS AND SERVICES PROVIDED FROM OUTSIDE**

Dry ice transport. Supplier selection and service assignment are carried out in accordance with applicable regulations.

7. **ASSIGNED VALUE, ROBUST RELATIVE STANDARD DEVIATION AND PERFORMANCE ASSESSMENT**

In order to minimize the influence of anomalous results on the statistical evaluation, the assigned value, or true value, is the average value calculated using robust statistics, as described in Annex C of ISO13528:2022 through algorithm A.

As established in point D.1.3.1 of Annex D of ISO 13528:2022, to proceed with the statistical evaluation a number of measurements greater than 12 is required ($n > 12$). If a sufficient number of data is not available, the statistical evaluation of the performance will not be carried out, the final report will report the data as is, without the z-score values, and adequate information will be provided that led to these conclusions.

	<p>Linea guida Fitofarmaci e Nitrati Pesticides and Nitrates Guideline</p>	<p>Allegato 6 I71201/PT rev.0</p>
---	--	---------------------------------------

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PT Pesticides

Gross errors mean errors of one or more orders of magnitude compared to the theoretical increment.

PT Nitrates

Gross errors mean errors of one or more orders of magnitude compared to the average value calculated when evaluating homogeneity.

If homogeneity is not experimentally assessed, the search for any gross errors will take place by comparing the data with respect to the median value of all the participants' results.

Any gross errors will be highlighted in the final report/data statistics tables, indicating the participant(s) who returned such data; furthermore, the exclusion from the calculation of the assigned value will be highlighted.

To correctly calculate the assigned value, the following chronological steps are performed:

1. evaluate, and eliminate, any grossly anomalous data,
2. evaluate that the number of measurements is greater than 12 ($n > 12$),
3. evaluate the possible presence of anomalous data (outliers) through the Minitab application or other statistical tests (Dixon or Grubbs test),
4. calculate the assigned value.

As established in point 6.5.2 of the ISO 13528:2022 Standard, algorithm A is able to compensate for the presence of anomalous data up to a percentage of no more than 20% of the entire data population. In the case of a percentage greater than 20%, the participants' performances will be calculated by applying parametric statistics, using the average of the data received from the participants as the assigned value, after eliminating anomalous data, and 25% of the average itself as the standard deviation. The z-scores will also be calculated for the eliminated outlier values and will be highlighted in the final report.

To evaluate the participants' performance, the z-score parameter is calculated through the mathematical relationship:

$$z = \frac{X_i - X_{(PT)}}{\sigma_{(PT)}}$$

Where:

- $x_{(i)}$: represents the value found by the laboratory for an active substance
- $x_{(pt)}$: represents the value found by the robust statistical processing for the PT
- $\sigma_{(pt)}$: represents the standard deviation FFP equal to 25% of the value assigned for PT Pesticides, and the FFP equal to 12% for PT Nitrates.

	<p>Linea guida Fitofarmaci e Nitrati Pesticides and Nitrates Guideline</p>	<p>Allegato 6 I71201/PT rev.0</p>
---	--	---------------------------------------

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The uncertainty associated with each parameter is calculated with the following equation:

$$u_{(XPT)} = 1.25 \times \frac{S^*}{\sqrt{n}}$$

Where:

- $u_{(XPT)}$: represents the uncertainty associated with the value assigned for each analyte.
- S^* : represents the robust standard deviation, calculated according to what is reported in Annex C of the ISO 13528:2022 standard.
- n : represents the number of measures.

a) The following relation is evaluated:

$$u_{(XPT)} < 0.3\sigma_{PT}$$

Where:

- $\sigma_{(PT)}$: standard deviation used in the evaluation of the PT equal to 25% of the value assigned for PT Plant Protection Products, and 12% of the value assigned for PT Nitrates.

b) The homogeneity of the objects is evaluated before shipping with the following relation:

$$S_s \leq 0.3\sigma_{om}$$

Where:

- S_s = standard deviation between samples, calculated as described in point B.3 of the ISO 13528:2022.
- σ_{om} = where FFP-RSD = 0.15 x average value of the homogeneity of each active substance present in the PT
- Control value: $0.3 \times \sigma_{om}$
- 0.3 = see Note 1 paragraph B.2.2 of ISO 13528:2022

c) The stability of the objects is evaluated with the following relation:

$$|x_i - y_i| \leq 0.3\sigma_{PT}$$

Where:

- x_i = average value of the samples on Day 1.
- y_i = mean value of samples on Day 2 or mean value of samples on Day 3.
- σ_{PT} = standard deviation used in the evaluation of the PT equal to 25% of the value assigned for the PT Pesticides, equal to 12% of the value assigned for the PT Nitrates
- 0.3: value obtained from Annex E.2 of ISO 13528:2022.

If one or more relation in points a); b); or c) are not respected, the z'-score will be calculated with the following relation:

$$z' = \frac{X_i - X_{PT}}{\sqrt{\sigma_{PT}^2 + u_{(XPT)}^2}}$$

Where:

- $x_{(i)}$: represents the value found by the laboratory for an active substance
- $x_{(pt)}$: represents the value found by the robust statistical processing for the PT
- $\sigma_{(pt)}$: represents the FFP standard deviation of the PT of interest
- $u_{(xpt)}$: the uncertainty associated with the assigned value for each analyte

The z'-score will be calculated even if homogeneity and/or stability is not experimentally evaluated.

	<p>Linea guida Fitofarmaci e Nitrati Pesticides and Nitrates Guideline</p>	<p>Allegato 6 I71201/PT rev.0</p>
---	--	---------------------------------------

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The meaning to be attributed to “z - score” and “z' – score” is interpreted according to the following table :

value of z	assessment
$ z \leq 2$	satisfying
$2 < z < 3$	questionable
$ z \geq 3$	unsatisfactory

8. STEPS TO FOLLOW

For the current year, participation in the proficiency tests involves 4 essential steps:

1. Each laboratory must complete and send to Arpae Annex 3 - participation form, received from the Provider via e-mail, by the date indicated in the communication. It must also clearly indicate the sessions and tests of interest. The e-mail address that is communicated through the participation form is the only one to be used for communications. Any changes must be formalized.
2. Participating laboratories must respect the times indicated for sending results for each session, using exclusively Google Sheets.
3. ARPAE evaluates the results once the deadline for submission has passed and **publishes a report for the pesticide test and a report for the nitrate test in PDF format on its website, on the proficiency testing page. A communication is sent to session participants, including a link to view the report.** This final report includes information on test preparation, homogeneity, stability, statistical evaluation and a graphic representation of the results of the participating laboratories, as well as any other relevant information.
4. Participating laboratories will receive the invoice related to costs for the current year by December 31, 2026, including costs for samples and shipping costs. The participating laboratories must pay what is due and reported. The costs are those reported in Annex 2, with 2026 costs.

9. SAMPLE QUANTITY

Participating laboratories receive at each session of interest:

- approximately 150 g of augmented sample for the determination of pesticides
- approximately 150 g of "blank" sample for the determination of pesticides
- approximately 50 g of sample for the determination of nitrates

	<p>Linea guida Fitofarmaci e Nitrati Pesticides and Nitrates Guideline</p>	<p>Allegato 6 I71201/PT rev.0</p>
---	--	---------------------------------------

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10. **SAMPLE SHIPPING**

All test samples are frozen and shipped in polystyrene boxes containing dry ice, which are then contained in cardboard boxes.

The expected shipping dates are those indicated in Annex 1, PT program for the year 2026 .

In the event of damage or loss of objects, the participating laboratory may request a second shipment without payment.

11. **WARNINGS ON SPECIMEN HANDLING**

Once received, test samples must be kept frozen until analysis, to avoid any deterioration. Once thawed, the samples must be shaken before taking the necessary portion.

All participants can use their usual analytical procedures and reference materials for identification and quantitation.

The sample should be handled as a normal routine sample.

Care must be taken when handling any dry ice still present in the package.

12. **RESULTS AND SIGNIFICANT FIGURES**

Results are returned via Google Sheet shared with each participant for each individual session. To view the parameters it is necessary to fill in the status of the sample upon arrival and the methods used for the determinations by selecting one of the items from the menu.

From the moment the sharing is opened until the deadline, it is possible to make an unlimited number of changes and/or additions, the Google Sheet automatically saves every change, so, there is no need to save the document.

The results must be expressed by typing the point as a decimal separator, it is not necessary to use any separator to identify the thousands.

No results should be reported if an active substance is not tested for, is not found upon analysis, or is found below the laboratory's LOQ. In these cases it is indicated as ND in the final paper.

No result should be reported if an active substance is also found in the "blank" sample. An active substance present in the test, analyzed and not detected, which corresponds to a participant's Limit of Quantification (LOQ) lower than the assigned value, is considered NR (Not Detected) and corresponds to a z-score of 5.

An active substance not present in the test, but detected, leads to an error for which the laboratory will be given a z-score of 5.

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In the performance evaluation, any result below the LOQ of the Proficiency test Provider (PTP), equal to 0.010 mg/kg, will not be considered.

In cases where the assigned value is less than a factor of 3 times the LOQ of the PTP, generally no false negatives or false positives will be assigned. The PTP may decide to take case-by-case decisions in this regard after considering all relevant factors such as the distribution of results and the LOQs of the laboratories concerned.

According to the DG-SANTE guidelines, an analyte result must be corrected for method bias if it exceeds 20%. Unless the method used inherently accounts for bias, such as:

- the use of stable isotope-labeled analogues of the target analytes as internal standards (ILIS), added to the analytical portion at an early stage in the procedure.**
- a "procedural calibration" approach.**
- a "standard-spiking" approach, where analytes are added to the analytical portions prior to extraction.**

Laboratories are required to report the result, expressed in mg/kg, with any recovery factor already applied. The percent recovery factor should not be reported in any cell of the Google Sheet used for data output.

Results with bias greater than 20% (apparent recovery outside the range of 80-120%), for which no bias correction was applied, may be omitted from the population used to calculate the assigned value.

The value of the concentrations of the active substances found during the analysis is indicated with a single figure in mg/kg and the number of significant figures is that indicated in the SANTE document in force, a significant figure may be added for statistical purposes.

conc interval . (mg/kg)	significant figures	examples
LOQ ≤ x < 10	2	0.053; 0.53; 5.3
x ≥ 10	3 or integer	35.2; 352

At the time of the deadline, sharing will be removed from the Google Sheet and the PDF file of their Google Sheet with the results entered will be sent to each individual participant.

13. APPEALS and COMPLAINTS

In accordance with the Arpae QMS, participants in the PTs have the possibility to present appeals or complaints to the Provider, by sending written communication via email to the email address proficiencytesting@arpae.it; using the form on the page dedicated to Proficiency tests on the Arpae website.

The nature of the complaint depends on how the service was provided: delays in the scheduled times, thawed or damaged samples.

	<p>Linea guida Fitofarmaci e Nitrati Pesticides and Nitrates Guideline</p>	<p>Allegato 6 I71201/PT rev.0</p>
---	--	---------------------------------------

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The appeal can be submitted after the final report is issued. The participant can contest the evaluation of their performance by presenting objective documentation of their motivations. Once the appeal or complaint has been received, Arpae communicates receipt to the participant and has 30 calendar days to respond via email and possibly amend the final report.

14. RESIDUAL OBJECTS

At the conclusion of each PT scheme, the excess objects are stored at a controlled freezing temperature as reported in the Arpae QMS, for a period of 30 calendar days from the date of issue of revision 0 of the final report. The objects remain available to the participants who have the right to request an additional amount if they wish to use the material for their own purposes.

It is specified that the stability of the material is guaranteed until the deadline for transmitting the results communicated in the email that is sent upon delivery of the samples. The participants are responsible for the logistics and cost of objects shipping.

15. STORAGE

All records relating to the PT schemes (e.g. results transmitted by participants, statistical calculations, reports, documents and any relevant communication) are conserved on computer by the Coordinator for at least 5 years. The support for computer storage is managed in accordance with P71301/ER and other specific instructions (backup).